

Package: photobiologyInOut (via r-universe)

October 24, 2026

Type Package

Title Read Spectral and Logged Data from Foreign Files

Version 0.4.33.9000

Date 2026-05-09

Description Functions for reading, and in some cases writing, foreign files containing spectral data from spectrometers and their associated software, output from daylight simulation models in common use, and some spectral data repositories. As well as functions for exchange of spectral data with other R packages. Part of the 'r4photobiology' suite, Aphalo P. J. (2015) [doi:10.19232/uv4pb.2015.1.14](https://doi.org/10.19232/uv4pb.2015.1.14).

License GPL (>= 2)

VignetteBuilder knitr

Depends R (>= 4.2.0), photobiology (>= 0.14.0)

Imports methods, tools, utils, stats, SunCalcMeeus (>= 0.1.3), stringr (>= 1.5.2), lubridate (>= 1.9.4), anytime (>= 0.3.12), tibble (>= 3.3.0), dplyr (>= 1.1.4), tidyr (>= 1.3.1), tidyselect (>= 1.2.1), readr (>= 2.1.5), readxl (>= 1.4.5), colorSpec (>= 1.8-0), jsonlite (>= 2.0.0)

Suggests spacesXYZ (>= 1.6-0), hyperSpec (>= 0.100.2), pavo (>= 2.9.0), fda.usc (>= 2.2.0), knitr (>= 1.50), rmarkdown (>= 2.29), ggplot2 (>= 4.0.0), ggspectra (>= 0.3.17), photobiologyWavebands (>= 0.5.4), testthat (>= 3.2.3)

LazyLoad yes

ByteCompile true

Encoding UTF-8

URL <https://docs.r4photobiology.info/photobiologyInOut/>

BugReports <https://github.com/aphalo/photobiologyinout/issues/>

Config/roxygen2/version 8.0.0

Config/pak/sysreqs libicu-dev libx11-dev

Repository <https://aphalo.r-universe.dev>
Date/Publication 2026-05-12 13:30:51 UTC
RemoteUrl <https://github.com/aphalo/photobiologyInOut>
RemoteRef HEAD
RemoteSha cb8f8c15eb14892fd10cd3f15b375908c7191a4a

Contents

photobiologyInOut-package	3
colorSpec2mspct	4
hyperSpec2mspct	6
mspct2colorSpec	7
mspct2fddata	9
qtuv_clouds	10
qtuv_s.e.irrad	11
read_asd_txt	14
read_ASTER_txt	17
read_avaspec_csv	19
read_cid_spectravue_csv	20
read_CIE_csv	23
read_csi_dat	24
read_fmi_cum	26
read_fmi2mspct	27
read_foreign2mspct	29
read_FReD_csv	30
read_li180_txt	31
read_licor_prn	33
read_macam_dta	35
read_oo_jazirrad	36
read_oo_ovirrad	38
read_oo_pidata	40
read_oo_ssirrad	42
read_qtuv_txt	44
read_spectrapen_csv	45
read_tuv_usrout	46
read_uvspec_disort	47
read_uvspec_disort_vesa	48
read_walz_lsa_xlsx	49
read_wasatch_csv	52
read_yoctopuce_csv	54
rspec2mspct	56
spct_CCT	57
spct_CRI	58
spct_SSI	59

Index

photobiologyInOut-package

photobiologyInOut: Read Spectral and Logged Data from Foreign Files

Description

Functions for reading, and in some cases writing, foreign files containing spectral data from spectrometers and their associated software, output from daylight simulation models in common use, and some spectral data repositories. As well as functions for exchange of spectral data with other R packages. Part of the 'r4photobiology' suite, Aphalo P. J. (2015) [doi:10.19232/uv4pb.2015.1.14](https://doi.org/10.19232/uv4pb.2015.1.14).

Data acquisition

The support for Ocean Insight, formerly Ocean Optics, spectrometers in package 'photobiologyInOut' is limited to the import of data acquired with Ocean Optics' software as is. In contrast, package 'ooacquire', part of these same suite, makes it possible to control, modify settings and acquire spectral data from Ocean Optics spectrometers directly from within R. 'ooacquire' also supports the conversion of raw-counts data into physical quantities.

Warning!

Most of the file formats supported are not standardized, and are a moving target because of changes in instrument firmware and support software. In addition the output format, especially with models, can depend on settings that users can alter. So do check that import is working as expected, and if not, please please raise an issue and upload one example of an incorrectly decoded file.

Note

From version 0.4.4 the time zone (tz) used for decoding dates and times in files imported defaults to "UTC". In most cases you will need to pass the tz (or the locale) where the file was created as an argument to the functions!

Author(s)

Maintainer: Pedro J. Aphalo <pedro.aphalo@helsinki.fi> ([ORCID](#))

Authors:

- Pedro J. Aphalo <pedro.aphalo@helsinki.fi> ([ORCID](#))

Other contributors:

- Titta K. Kotilainen ([ORCID](#)) [contributor]
- Glenn Davis <gdavis@gluonics.com> [contributor]

References

Aphalo, Pedro J. (2015) The r4photobiology suite. UV4Plants Bulletin, 2015:1, 21-29. [doi:10.19232/uv4pb.2015.1.14](https://doi.org/10.19232/uv4pb.2015.1.14).

See Also

Useful links:

- <https://docs.r4photobiology.info/photobiologyInOut/>
- Report bugs at <https://github.com/aphalo/photobiologyinout/issues/>

colorSpec2mspct	<i>Convert 'colorSpec::colorSpec' objects</i>
-----------------	---

Description

Convert 'colorSpec::colorSpec' objects into spectral objects (xxxx_spct, xxxx_mspct) as defined in package 'photobiology' and vice versa preserving as much information as possible.

Usage

```
colorSpec2mspct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.source_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.source_mspct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.response_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.response_mspct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.filter_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.filter_mspct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.reflector_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.reflector_mspct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.chroma_mspct(x, multiplier = 1, ...)

colorSpec2spct(x, multiplier = 1, ...)
```

```

colorSpec2chroma_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.chroma_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.chroma_mspct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.generic_spct(x, multiplier = 1, ...)

## S3 method for class 'colorSpec'
as.generic_mspct(x, multiplier = 1, ...)

```

Arguments

<code>x</code>	colorSpec object
<code>multiplier</code>	numeric A multiplier to be applied to the 'spc' data to do unit or scale conversion.
<code>...</code>	currently ignored.

Details

Objects of class `colorSpec::colorSpec` do not contain metadata or class data from which the units of expression could be obtained. When using function `colorSpec2mspct` the user needs to use parameter `multiplier` to convert the data to what is expected by the object constructors defined in package 'photobiology' but should only rarely need to use parameter `spct.data.var` to select the quantity.

`colorSpec::colorSpec` objects may use memory more efficiently than spectral objects of the classes for collections of spectra defined in package 'photobiology' as wavelengths are assumed to be the same for all member spectra, and stored only once while this assumption is not made for collections of spectra, allowing different wavelengths and lengths for the component spectra. When using `as.colorSpec` methods to convert collections of spectra into `colorSpec` objects, if the wavelengths of the individual spectra differ, only the shared range of wavelengths is retained and within this range, wavelength values are made consistent by interpolation.

Warning!

Always check the sanity of the imported or exported data values, as guessing is needed when matching the different classes, and the functions defined here are NOT guaranteed to return valid data without help from the user through optional function arguments.

Note

In `generic_mspct` objects, wavelengths are stored for each spectrum, individual `generic_spct` objects. However, as spectral classes are derived from 'tbl_df' in many cases no redundant copies of wavelength data will be made in memory in spite of the more flexible semantics of the objects.

Examples

```
# example run only if 'colorSpec' is available
if (requireNamespace("colorSpec", quietly = TRUE)) {
  library(colorSpec)
  colorSpec2mspct(Fs.5nm)
  colorSpec2spct(Fs.5nm)
  colorSpec2mspct(C.5nm)
  colorSpec2spct(C.5nm)
}
```

hyperSpec2mspct

Convert 'hyperSpec::hyperSpec' objects

Description

Convert `hyperSpec::hyperSpec` objects containing VIS and UV radiation data into spectral objects (`xxxx_spct`, `xxxx_mspct`) as defined in package 'photobiology' and vice versa, preserving as much information as possible. As `hyperSpec` can contain other kinds of spectral data, it does make sense to use these functions only with objects containing data that can be handled by both packages.

Usage

```
hyperSpec2mspct(x, member.class, spct.data.var, multiplier = 1, ...)
```

```
hyperSpec2spct(x, multiplier = 1, ...)
```

```
mspct2hyperSpec(x, spct.data.var, multiplier = 1, ...)
```

```
spct2hyperSpec(x, spct.data.var = NULL, multiplier = 1, ...)
```

Arguments

<code>x</code>	hyperSpec object
<code>member.class</code>	character One of the spectrum classes defined in package 'photobiology'.
<code>spct.data.var</code>	character The name to be used for the 'spc' data when constructing the spectral objects.
<code>multiplier</code>	numeric A multiplier to be applied to the 'spc' data to do unit or scale conversion. For example "a.u." units in some examples in package 'hyperSpec' seem to have scale factors applied.
<code>...</code>	currently ignored.

Warning!

Always check the sanity of the imported or exported data values, as guessing is needed when matching the different classes, and the functions defined here are NOT guaranteed to return valid data without help from the user through optional function arguments.

Note

Objects of class `hyperSpec::hyperSpec` contain metadata or class data from which the quantity measured and the units of expression can be obtained. However, units as included in the objects are not well documented making automatic conversion difficult. When using this function the user may need to use parameter `multiplier` to scale the data to what is expected by the object constructors defined in package 'photobiology' and use parameter `spct.data.var` to select the quantity.

`hyperSpec::hyperSpec` objects may use memory more efficiently than spectral objects of the classes for collections of spectra defined in package 'photobiology' as wavelengths are assumed to be the same for all member spectra, and stored only once while this assumption is not made for collections of spectra, allowing different wavelengths and lengths for the component spectra. Wavelengths are stored for each spectrum, but as spectral classes are derived from 'tbl_df' in many cases no redundant copies of wavelength data will be made in memory in spite of the more flexible semantics of the objects.

Examples

```
# example run only if 'hyperSpec' is available
if (requireNamespace("hyperSpec", quietly = TRUE)) {
  library(hyperSpec)
  data(laser)
  wl(laser) <-
    list(wl = 1e7 / (1/405e-7 - wl(laser)),
         label = expression (lambda / nm))
  laser.mspct <- hyperSpec2mspct(laser, "source_spct", "s.e.irrad")
  class(laser.mspct)
}
```

mspct2colorSpec	<i>Convert into 'colorSpec::colorSpec' objects</i>
-----------------	--

Description

Convert spectral objects (`xxxx_spct`, `xxxx_mspct`) as defined in package 'photobiology' into color-Spec objects preserving as much information as possible.

Usage

```
mspct2colorSpec(x, spct.data.var = NULL, multiplier = 1, ...)

spct2colorSpec(x, spct.data.var = NULL, multiplier = 1, ...)

chroma_spct2colorSpec(x, spct.data.var = NULL, multiplier = 1, ...)

## S3 method for class 'generic_mspct'
as.colorSpec(x, spct.data.var = NULL, multiplier = 1, ...)
```

```
## S3 method for class 'generic_spct'
as.colorSpec(x, spct.data.var = NULL, multiplier = 1, ...)

## S3 method for class 'chroma_spct'
as.colorSpec(x, spct.data.var = NULL, multiplier = 1, ...)
```

Arguments

<code>x</code>	R object
<code>spct.data.var</code>	character The name of the variable to read spectral data from.
<code>multiplier</code>	numeric A multiplier to be applied to the 'spc' data to do unit or scale conversion.
<code>...</code>	currently ignored.

Warning!

Always check the sanity of the returned data values, as guessing is needed when matching the different classes, and the functions defined here are NOT guaranteed to return valid data without help from the user through optional function arguments.

Note

Objects of class `colorSpec::colorSpec` do not contain metadata or class data from which the units of expression could be obtained. When using this function the user needs to use parameter `multiplier` to convert the data to what is expected by the object constructors defined in package 'photobiology' but should only rarely need to use parameter `spct.data.var` to select the quantity.

`colorSpec::colorSpec` objects may use memory more efficiently than spectral objects of the classes for collections of spectra defined in package 'photobiology' as wavelengths are assumed to be the same for all member spectra, and stored only once while this assumption is not made for collections of spectra, allowing different wavelengths and lengths for the component spectra. Wavelengths are stored for each spectrum, but as spectral classes are derived from 'tbl_df' in many cases no redundant copies of wavelength data will be made in memory in spite of the more flexible semantics of the objects.

Examples

```
if (requireNamespace("colorSpec", quietly = TRUE)) {
  library(colorSpec)
  as.colorSpec(polyester.spct)
  as.colorSpec(sun.spct)
  as.colorSpec(filter_mspct(list(polyester.spct, yellow_gel.spct)))
}
```

mspct2fdata

*Convert spectra into 'fda.usc::fdata' objects***Description**

Convert spectral objects (xxxx_spct, xxxx_mspct) as defined in package 'photobiology' into `fda.usc::fdata` objects, preserving as much information as possible. As `fdata` objects can contain other kinds of data, the reverse conversion is supported (experimentally) and mainly for 'fdata' objects returned by the functional data analysis methods from package `fda.usc` to spectral data previously exported in the opposite direction.

Usage

```
mspct2fdata(x, spct.data.var = NULL, multiplier = 1, ...)
```

```
spct2fdata(x, spct.data.var = NULL, multiplier = 1, ...)
```

```
fdata2spct(x, multiplier = 1, member.class = NULL, drop.idx = FALSE, ...)
```

```
fdata2mspct(x, multiplier = 1, member.class = NULL, drop.idx = FALSE, ...)
```

Arguments

<code>x</code>	generic_mspct or generic_spct object or an object belonging to a derived class, or an object of class 'fdata' depending on the function.
<code>spct.data.var</code>	character The name of the column containing data to export. If NULL the first spectral data column found is used.
<code>multiplier</code>	numeric A multiplier to be applied to the 'spc' data to do unit or scale conversion.
<code>...</code>	possibly additional named arguments passed to object constructors.
<code>member.class</code>	character Name of the class of the spectrum or of the members of the collection of spectra.
<code>drop.idx</code>	logical Flag indicating whether to drop or keep <code>idx.var</code> in the collection members.

Warning!

When converting multiple spectra, all the spectra to be included in the `fdata` object must share the same wavelength values. Spectra that do not fulfil this condition will be skipped. The data variable needs also to be present in all individual spectra as no conversions are applied automatically by this function. If a different name, indicating a different quantity or a different base of expression is encountered, the affected spectrum is skipped with a warning.

Examples

```

if (requireNamespace("fda.usc", quietly = TRUE)) {
  # from spectra to fdata
  sun.fdata <- spct2fdata(sun.spct)
  str(sun.fdata)
  polyester.fdata <- spct2fdata(polyester.spct)
  str(polyester.fdata)
  # from fdata to spectra
  fdata2spct(sun.fdata)
  fdata2spct(sun.fdata, drop.idx = TRUE)
  fdata2spct(polyester.fdata, drop.idx = TRUE)
}

```

qtuv_clouds

Clouds descriptor

Description

Constructor of a named list of parameter values to be used as argument to parameter clouds of function `qtuv_s.e.irrad()`.

Usage

```
qtuv_clouds(cloud.type = "clear.sky")
```

Arguments

`cloud.type` character One of "clear.sky", "cirrus", "stratocumulus" or "overcast".

Details

This function provide a rough approximation for parameter values. In reality there is large variation in the cloud optical depths (COD) and in the elevation at which clouds are located, within each type of cloud. The TUV model assumes a continuous uniform cloud layer, thus the normally discontinuous cover of cumulus clouds cannot be described.

Value

A one-row data frame with members "optical.depth", "base", "top" and "label".

Examples

```

qtuv_clouds("clear.sky")
qtuv_clouds("cirrus")
qtuv_clouds(c("clear.sky", "cirrus"))

```

qtuv_s.e.irrad

Spectral irradiance from the Quick TUV calculator

Description

Call the Quick TUV calculator web server and return a `source_spct` object with the simulated spectral energy irradiance data.

Usage

```
qtuv_s.e.irrad(
  w.length = list(wStart = 280, wStop = 420, wIntervals = 140),
  sun.elevation = NULL,
  geocode = data.frame(lon = 0, lat = 51.5, address = "Greenwich"),
  time = lubridate::now(),
  tz = NULL,
  locale = readr::default_locale(),
  ozone.du = 300,
  albedo = 0.1,
  ground.altitude = 0,
  measurement.altitude = NULL,
  clouds = data.frame(optical.depth = 0, base = 4, top = 5),
  aerosols = data.frame(optical.depth = 0.235, ssaaer = 0.99, alpha = 1),
  num.streams = 2,
  spectra = list(direct = 1, diffuse.down = 1, diffuse.up = 0),
  added.vars = NULL,
  label = "",
  server.url = "https://www.acom.ucar.edu/cgi-bin/acom/TUV/V5.3/tuv",
  file = NULL
)

qtuv_m_s.e.irrad(
  w.length = list(wStart = 280, wStop = 420, wIntervals = 140),
  sun.elevation = NULL,
  geocode = data.frame(lon = 0, lat = 51.5, address = "Greenwich"),
  time = lubridate::now(),
  tz = NULL,
  locale = readr::default_locale(),
  ozone.du = 300,
  albedo = 0.1,
  ground.altitude = 0,
  measurement.altitude = NULL,
  clouds = data.frame(optical.depth = 0, base = 4, top = 5),
  aerosols = data.frame(optical.depth = 0.235, ssaaer = 0.99, alpha = 1),
  num.streams = 2,
  spectra = list(direct = 1, diffuse.down = 1, diffuse.up = 0),
  added.vars = NULL,
```

```

    label = "",
    server.url = "https://www.acom.ucar.edu/cgi-bin/acom/TUV/V5.3/tuv",
    file = NULL
  )

```

Arguments

<code>w.length</code>	list of parameters describing the wavelengths, or a numeric vector from which the parameters will be constructed.
<code>sun.elevation</code>	numeric Angle in degrees above the horizon. If NULL its value is computed from geocode and time, otherwise arguments passed to these two parameters are ignored.
<code>geocode</code>	data.frame with variables lon and lat as numeric values (degrees), and character variable address; nrow > 1, allowed for collections.
<code>time</code>	A "vector" of POSIXct time, with any valid time zone (TZ) is allowed, default is current time.
<code>tz</code>	character Time zone is by default read from the file.
<code>locale</code>	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
<code>ozone.du</code>	numeric Ozone column in Dobson units.
<code>albedo</code>	numeric Surface albedo (= reflectance) as a fraction of one.
<code>measurement.altitude, ground.altitude</code>	numeric Altitudes above sea level expressed in km.
<code>clouds</code>	data.frame Parameters optical.depth (vertical), top and base expressed in km; nrow > 1, allowed for collections.
<code>aerosols</code>	data.frame Parameters optical.depth (total extinction), ssaaer (cloud single scattering albedo) and alpha (wavelength dependence of optical depth); nrow > 1, allowed for collections.
<code>num.streams</code>	integer Number of streams used in computations, 2 or 4.
<code>spectra</code>	named list with weights for the different components of the spectrum.
<code>added.vars</code>	character vector Accepted member values are "sun.elevation", "zenith.angle", "date" and "ozone.du".
<code>label</code>	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
<code>server.url</code>	character The URL used to access the Quick TUV calculator server.
<code>file</code>	character The name under which the file returned by the server is locally saved. If NULL a temporary file is used and discarded immediately. File paths are supported when valid.

Details

The Quick TUV calculator, is an on-line freely accessible server running the TUV atmospheric chemistry and radiation transfer model with a simplified user interface. In this case, version 5.3 is called passing the parameter values passed as arguments in the call to `qtuv_s.e.irrad()`. The response is saved in a text file that is subsequently read with function `read_qtuv_txt()` into a `source_spct` object.

Function `qtuv_m_s.e.irrad()` calls `qtuv_s.e.irrad()` repeatedly accepting a numeric vector longer than one as argument, for parameters: `sun.elevation`, `time` or `ozone.du`, and data frames with `nrow >= 1`. In a given call, only one parameter at a time can obey multiple values, with others currently truncated to the first value.

The formal parameter names are informative and consistent with other functions in the R for Photobiology Suite and differ from the short names used for the parameters in the FORTRAN code of the TUV model. In the case of `w.length` two ways of specifying wavelengths are supported. Some defaults also differ from those of the Quick TUV calculator.

In the current implementation, `qtuv_m_s.e.irrad()`, accepts multiple values as arguments for only one parameter at a time. In the case of elevation, both `'ground.elevation'` and `'measurement.elevation'` can have each one or more values. When too many values are passed in the call, only the first one is used.

Value

In the case of `qtuv_s.e.irrad()`, a `source_spct` object obtained by finding the center of each wavelength interval in the Quick TUV output file, and adding the variables listed in `added.vars`. In the case of `qtuv_m_s.e.irrad()`, a `source_mspct` object containing a collection of such spectra.

Side effect

If a file name is passed as argument, the data as downloaded are saved into persistent files, one file per spectrum. The names of the saved files always end in `'txt'`.

Warning!

This function connects to a server managed by UCAR, the University Corporation for Atmospheric Research located in the U.S.A. to obtain simulated spectral data. UCAR manages the U.S. National Science Foundation National Center for Atmospheric Research (NSF NCAR) on behalf of NSF. As any download with the HTTPS protocol, using this function entails some risk. To minimize the risk, the returned page is saved as plain text, checked for conformity with the expected content, and if valid decoded into an R data object. When using the default argument `file = NULL`, the file used is a temporary one and is deleted before the function returns the call, irrespective of it being conformant or not.

The administrators of the Quick TUV Calculator at UCAR suggest a maximum load of approximately 100 spectral simulations per day and user. For larger workloads they encourage the local installation of the TUV model which is open-source and freely available. A local installation, also allows access to the full set of input parameters and outputs. Currently a local instance of the TUV model can be called from R with package `'foqat'`.

Note

The Quick TUV calculator has multiple output modes that return different types of computed values. The use of output mode 5 is hard-coded in this function as other modes return summary values rather than spectral data. Package 'foqat' provides a more flexible alternative supporting other output modes in addition to mode 5.

If the argument passed to `w.length` is a numeric vector, the range and length are used to reconstruct the accepted parameters. The returned spectrum has always a uniformly spaced wavelengths.

When using this function, more detailed metadata are available than when reading an output file, as not all the simulation input parameters are listed in the output text.

In interactive use of the Quick TUV Calculator, the same parameters as accepted by `qtuv$.e.irradiance()` as arguments are entered through the web interface at https://www.acom.ucar.edu/Models/TUV/Interactive_TUV/. This page together with its documentation, can be consulted for additional information on the parameters and the model.

References

https://www.acom.ucar.edu/Models/TUV/Interactive_TUV/. This URL could change in the future as well as the server URL. The formal parameter `server.url` was included only for use in such a case.

Sasha Madronich (2017-2021) Tropospheric Ultraviolet and Visible radiation (TUV) model. <https://www2.acom.ucar.edu/modeling/tropospheric-ultraviolet-and-visible-tuv-radiation-model>. Visited on 2024-08-29.

read_asd_txt

Read File Saved by ASD's conversion tool

Description

Read the spectral data and parse the header of a energy irradiance, a reflectance, a transmittance data, a raw-detector-counts data or arbitrary data file as output by the text conversion tool for ASD spectrometers when set for tab-, comma- or semicolon-separated output. ASD's field spectrometers measure VIS and SWIR radiation.

Usage

```
read_asd_txt(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  s.qty = NULL,
  scale.factor = 1,
  range = NULL
)
```

Arguments

file	character string Path to the file to be read, following R's use of forward slashes as separator for folder names.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat, and optionally address used to set attribute "where.measured".
label	character string to which to set the "what.measured" attribute. If NULL the value of basename(file) is used, and if NA the "what.measured" attribute is not set.
tz	character A time zone recognized by R. If NULL it is extracted from locale. If "", the default, the local time zone is used.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), with time in UTC but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names. Its value must match that used to write the imported file, which is not necessarily the default one or the local one.
s.qty	character string The name of the quantity to be read as by object constructors from 'photobiology'. If NULL quantity is guessed from file header.
scale.factor	numeric Multiplier to be applied to the spectral data so that it matches s.qty in base of expression and units. In most cases the default of 1.0 is correct.
range	a numeric vector of length two, or any other object for which function range() will return a range of wavelengths expressed in nanometres.

Details

The header of the file is first decoded and parsed to extract the time of data acquisition and serial number of the spectrometer, to locate the start of the spectral data and to detect the column separator used. The time in the header is the local time with no time zone information. Thus with the default `tz = ""` the time is re-expressed as UTC with a correction. When file import takes place on a different time zone than the measurement the measurement time zone must be supplied by the user as an argument to parameter `tz` or as part of the `locale` after setting `tz = NULL`. The metadata fields in the header are located by text matching.

Spectral irradiance is returned as an object of class [source_spct](#), reflectance as [reflector_spct](#), transmittance as [filter_spct](#), raw detector counts as [raw_spct](#), and arbitrary spectral values as [generic_spct](#) objects with metadata stored in attributes `when.measured`, `what.measured`, and `how.measured`. The spectral quantity is guessed from the header metadata by default, but an argument to `s.qty` can be passed to override this default.

The value stored in the `how.measured` attribute includes the spectrometer serial number extracted from the file header. If an argument is passed to parameter `geocode`, its value is saved in attribute `where.measured` (currently the `geocode` is not extracted from the file header). The object's comment always gives a text that includes the file name, time of import, function name and the version of packages 'photobiology' and 'photobiologyInOut' used. The file header in whole is copied to attribute `file.header`.

Value

A `source_spct` object with a column `s.e.irrad` with spectral energy irradiance in $Wm^{-2}nm^{-1}$, or a `reflector_spct` object with a column `Rfr` with fractional reflectance, or a `filter_spct` object with a column `Tfr` with fractional total transmittance. In all cases a column `w.length` contains wavelengths in nanometres and the attributes `comment`, `what.measured`, `when.measured`, `how.measured`, `where.measured` and `file.header` containing metadata for the spectrum.

Note

Some files with ASD spectral data contain nul ASCII characters in the header. These are ignored with a warning while the data are still read. Depending on the program used to view these data files, only the top three lines will be displayed, or the file will be detected as containing binary data and not displayed.

References

<https://www.malvernpanalytical.com/en/>

Examples

```
# energy spectral irradiance file

file.name <-
  system.file("extdata", "asd-e-irrad-sky.tsv",
             package = "photobiologyInOut",
             mustWork = TRUE)

asd.source_spct <-
  read_asd_txt(file = file.name,
              tz = "Europe/Helsinki")

class_spct(asd.source_spct)
summary(asd.source_spct)
getWhenMeasured(asd.source_spct)
getWhatMeasured(asd.source_spct)
getHowMeasured(asd.source_spct)
cat(comment(asd.source_spct))
e_irrad(asd.source_spct)
q_irrad(asd.source_spct, w.band = c(400, 700),
       scale.factor = 1e6)

asd_clipped.source_spct <-
  read_asd_txt(file = file.name,
              tz = "Europe/Helsinki",
              range = c(400, 700))

class_spct(asd_clipped.source_spct)
summary(asd_clipped.source_spct)

# spectral reflectance file
```

```

file.name <-
  system.file("extdata", "asd-Rfr-panel-50pc.tsv",
             package = "photobiologyInOut",
             mustWork = TRUE)

asd.reflector_spct <-
  read_asd_txt(file = file.name,
             tz = "Europe/Helsinki")

class_spct(asd.reflector_spct)
summary(asd.reflector_spct)
getWhenMeasured(asd.reflector_spct)
getWhatMeasured(asd.reflector_spct)
getHowMeasured(asd.reflector_spct)
cat(comment(asd.reflector_spct))

asd_clipped.reflector_spct <-
  read_asd_txt(file = file.name,
             tz = "Europe/Helsinki",
             range = c(400, 700))

class_spct(asd_clipped.reflector_spct)
summary(asd_clipped.reflector_spct)

# Raw-counts data

file.name <-
  system.file("extdata", "asd-raw-gravel.tsv",
             package = "photobiologyInOut",
             mustWork = TRUE)

asd.raw_spct <-
  read_asd_txt(file = file.name,
             tz = "Europe/Helsinki")

class_spct(asd.raw_spct)
summary(asd.raw_spct)
getWhenMeasured(asd.raw_spct)
getWhatMeasured(asd.raw_spct)
getHowMeasured(asd.raw_spct)
cat(comment(asd.raw_spct))

```

read_ASTER_txt

Read File downloaded from ASTER data base.

Description

Reads and parses the header of a test file as available through the ASTER reflectance database. The Name field is retrieved and copied to attribute "what.measured". The header of the file is preserved as a comment.

Usage

```
read_ASTER_txt(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  npixels = Inf
)
```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Ignored.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
npixels	integer Number of pixels in spectral data.

Value

A raw_spect object.

Note

The header in these files has metadata information, but mostly on the origin of the data. For a date and/or geocode are to be added to the return object it must be supplied by the user. as well as the date-time. Some metadata is extracted and added as attributes, while the whole header is copied to the comment attribute.

References

<https://speclib.jpl.nasa.gov>

Baldrige, A.; Hook, S.; Grove, C. & Rivera, G. (2009) The ASTER spectral library version 2.0. Remote Sensing of Environment. 113, 711-715

Examples

```

file.name <-
  system.file("extdata", "aster-Rpc-drygrass.txt",
             package = "photobiologyInOut", mustWork = TRUE)

fred.spct <- read_ASTER_txt(file = file.name, npixels = Inf)

fred.spct
getWhatMeasured(fred.spct)
cat(comment(fred.spct))

```

read_avaspec_csv	<i>Read '.csv' File Saved by Avantes' Software for AvaSpec.</i>
------------------	---

Description

Reads and parses the header of a processed data file as output by the program Avaspec and then imports wavelength and spectral irradiance values. The file header has little useful metadata information.

Usage

```

read_avaspec_csv(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)

read_avaspec_xls(
  path,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)

```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, attribute when.measured is not set as the file header does not contain date and time information.

geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone used for interpreting times saved in the file header.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
path	Path to the xls/xlsx file

Value

A source_spct object.

References

<https://www.avantes.com/>

Examples

```
file.name <-
  system.file("extdata", "avantes-avaspec-e-irrad.csv",
             package = "photobiologyInOut", mustWork = TRUE)

avaspec.spct <- read_avaspec_csv(file = file.name)

avaspec.spct
what_measured(avaspec.spct)
when_measured(avaspec.spct)
how_measured(avaspec.spct)
cat(comment(avaspec.spct))
```

read_cid_spectravue_csv

Read File Saved by CID's SpectraVue.

Description

Read wavelength and spectral data from Measurements.CSV files exported from CID Bio-Sciences' **SpectraVue CI-710s** (not the older CI-710!) leaf spectrometer, importing them into R. Available metadata is also extracted from the file. Function `read_cid_spectravue_csv()` only accepts "row oriented" CSV files. These may contain multiple spectra, one per row.

Usage

```
read_cid_spectravue_csv(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  range = c(380, 1100),
  simplify = TRUE,
  absorbance.to = "object",
  strict.range = NA,
  ...
)
```

Arguments

file	character string Path to the file to be read, following R's use of forward slashes as separator for folder names.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat, and optionally address used to set attribute "where.measured".
label	character string to which to set the "what.measured" attribute. If NULL the value of basename(file) is used, and if NA the "what.measured" attribute is not set.
tz	character A time zone recognized by R. If NULL, the default, it is extracted from 'locale'.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names. Its value must match that used to write the imported file, which is not necessarily the default one or the local one.
range	a numeric vector of length two, or any other object for which function range() will return a range of wavelengths expressed in nanometres.
simplify	logical If TRUE, single spectra are returned as individual spectra instead of collections of length one.
absorbance.to	character Affects only absorbance measurements: "object", "all" or "A".
strict.range	logical Flag indicating whether off-range values result in an error (TRUE) instead of a warning (FALSE), or the test is disabled (NA).
...	additional arguments passed to the constructor of the 'filter_spect' object.

Details

SpectraVue's row-wise spectral Measurements.CSV files contain columns with metadata on the left edge, followed by columns with data for each of the 2048 pixels or wavelengths. The value in

column "Mode" indicates the quantity measured, decoded into variables Tpc, Rpc or A. The data in each row of the CSV file are read and stored in a `filter_spct`, `reflector_spct` or `object_spct` object. These objects are collected into a single `filter_mspct`, `reflector_mspct`, `object_mspct` or `generic_spct` object and returned.

Spectral data outside the range 400 nm to 1000 nm are very noisy and thus outside the valid range for the measurements. Out-of-range spectral data can also be caused by calibration drift. Consequently, reading of data is done always with the range check disabled, while whether a check is used before returning the collection of spectra depends on the argument passed to `strict.range` which by default is set to disable checks. This is done, because in most cases measurements from this instrument tend to require further processing before they comply with theoretical expectations of $Tfr + Rfr + Afr = 1$.

Value

An object of class `filter_spct`, `reflector_spct`, `object_spct` or `generic_mspct`.

Internal vs. total transmittance and absorbance

Spectravue returns transmittance values labelled with Transmittance as Mode. Transmittance values are not *total* as most of the scattered light transmitted is not detected. Absorbance (Abs) values returned labelled with Absorbance as Mode are for absorbance computed from the Transmittance. This estimate of absorbance overestimates real absorbance in the case of scattering materials like plant leaves. It is best to save spectral data acquired as absorbance into objects of class `'object_spct'` containing reflectance and transmittance, the default, as this preserves all the acquired data.

Specular vs. total reflectance

This function assumes that SpectraVue returns close to *total* reflectance readings. Given the optics of the instrument this is likely only a weak approximation.

Warning!!

CID's support has answered on 2022-05-19 that the extremely biased (plainly wrong!) values of transmittance measured by this instrument are due to a design flaw and that they are working on a solution for the problem. In practice, reflectance seems biased but usable as an instrument-specific quantity with arbitrary units. Transmittance and absorbance seem useless as values are wrong by about an order of magnitude. As of 2025-06-29 they are still to contact me again.

Note

SpectraVue creates three or four .CSV files for each series of measurements saved. Of these files, this function reads the one with name ending in `Measurements.CSV`. The first part of the file name gives the time of the session, but as the files can contain multiple spectra measured at different times, the time metadata is extracted separately for each spectrum. We provide a default argument for `range` that discards data for short and long wavelengths because values outside this range are according to the instrument's manual outside the usable range and in practice extremely noisy.

References

<https://cid-inc.com/>

Examples

```
# read file containing a single reflectance spectrum

file.name <-
  system.file("extdata", "cid-spectravue-Rpc-Measurements.csv",
             package = "photobiologyInOut", mustWork = TRUE)

cid_filter.spct <-
  read_cid_spectravue_csv(file = file.name)
summary(cid_filter.spct)

cid_filter.spct <-
  read_cid_spectravue_csv(file = file.name, simplify = FALSE)
summary(cid_filter.spct)

# read data measured as absorbance (A, Rpc and Tpc)

file.name <-
  system.file("extdata", "cid-spectravue-Abs-Measurements.csv",
             package = "photobiologyInOut", mustWork = TRUE)
cid.object_spct <-
  read_cid_spectravue_csv(file = file.name)
summary(cid.object_spct)

cid_A.filter_spct <-
  read_cid_spectravue_csv(file = file.name, absorbance.to = "A")
summary(cid_A.filter_spct)
```

read_CIE_csv

*Read '.CSV' files from CIE***Description**

Reads a CSV spectral data file and its companion JSON file with metadata as published by *International Commission on Illumination* (CIE) and then imports wavelengths and spectral values into one of the classes for spectral data defined in package 'photobiology'.

Usage

```
read_CIE_csv(file.name, label = NULL, simplify = FALSE)
```

Arguments

file.name	character string Name of CSV file.
label	character string. If NULL metadata read from the JSON file is used, and if NA the "what.measured" attribute is not set, and if a character string is passed, it is used to set the "what.measured" attribute.
simplify	logical If TRUE and the read file contained a single spectrum, extract the spectral object from the collection.

Details

The CSV file contains only numbers encoded as character strings, and the JSON file contains extensive metadata. The type of spectral data is encoded as part of the file name. If the original file name of the CSV file is passed as argument to parameter `file`, the function can retrieve all data and metadata, enough to return an R object of the correct class. The JSON file must be located in the same folder as the CSV file.

Value

Depending on the contents of the file, a `source_mspct` or `source_spct` object, a `response_spct` or `response_erspct` object, or a `chroma_spct` object, containing both data and metadata. The whole metadata JSON converted into an R 'list' is stored in attribute "header".

References

<https://cie.co.at/data-tables>

Examples

```
file.name <-
  system.file("extdata", "CIE_illum_C.csv",
             package = "photobiologyInOut", mustWork = TRUE)

CIE_illum_C.spct <- read_CIE_csv(file.name)
CIE_illum_C.spct

file.name <-
  system.file("extdata", "CIE_sle_photopic.csv",
             package = "photobiologyInOut", mustWork = TRUE)

CIE_sle_photopic.spct <- read_CIE_csv(file.name)
CIE_sle_photopic.spct
```

read_csi_dat

Read '.DAT' file(s) saved by modern Campbell Scientific loggers.

Description

Reads and parses the header of a processed data file as output by the PC400 or PC200W programmes extracting variable names, units and quantities from the header. Uses the comment attribute to store the metadata.

Usage

```
read_csi_dat(
  file,
  geocode = NULL,
  label = NULL,
  data_skip = 0,
  n_max = Inf,
  locale = readr::default_locale(),
  na = c("", "NA", "NAN"),
  ...
)
```

Arguments

file	Path to file as a character string.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
data_skip	integer Number of records (rows) to skip from the actual data block.
n_max	integer Maximum number of records to read.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
na	character Vector of strings to interpret as missing values. Set this option to character() to indicate no missing values.
...	Further named arguments currently passed to read_csv().

Value

read_csi_dat() returns a `tibble::tibble` object.

Note

This function is not useful for .DAT and .PRN files from old CSI loggers and software. Those were simple files, lacking metadata, which was stored in separate .FLD files.

References

<https://www.campbellsci.eu/>

read_fmi_cum	<i>Read daily cummulated solar spectrum data file(s).</i>
--------------	---

Description

Read one or more cumulated daily spectral irradiance file as output by Anders Lindors' model based on libRadTrans. Dates are read from the file header and parsed with the function supplied as `date.f`.

Usage

```
read_fmi_cum(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = "UTC",
  locale = readr::default_locale(),
  .skip = 3,
  .n_max = -1,
  .date.f = lubridate::ymd
)
```

```
read_m_fmi_cum(
  files,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = "UTC",
  .skip = 3,
  .n_max = -1,
  .date.f = lubridate::ymd
)
```

Arguments

<code>file</code>	Either a path to a file, a connection, or literal data (either a single string or a raw vector).
<code>date</code>	a <code>POSIXct</code> object to use to set the <code>"when.measured"</code> attribute. If <code>NULL</code> , the default, the date is extracted from the file header.
<code>geocode</code>	A data frame with columns <code>lon</code> and <code>lat</code> used to set attribute <code>"where.measured"</code> .
<code>label</code>	character string, but if <code>NULL</code> the value of <code>file</code> is used, and if <code>NA</code> the <code>"what.measured"</code> attribute is not set.
<code>tz</code>	character Time zone used for interpreting times saved in the file header.
<code>locale</code>	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that

	controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
.skip	Number of lines to skip before reading data—i.e. the number of rows in the header.
.n_max	Maximum number of records to read.
.date.f	A function for extracting a date-time from the file header passed as character string to its first argument and which returns a POSIXct object.
files	list or vector of paths each one with the same requirements as described for argument file.

Value

read_fmi_cum() returns a source_spct object with time.unit attribute set to "day" and when.measured attribute set to the date-time extracted from the header at the top of the read file.

read_m_fmi_cum returns a source_mspct containing one source_spct object for each one of the multiple files read.

Note

See [read_table](#) for details of acceptable values for file.

Examples

```
file.name <- system.file("extdata", "fmi-al-2014-08-21_cum.hel",
                        package = "photobiologyInOut", mustWork = TRUE)
fmi.spct <- read_fmi_cum(file = file.name)
```

read_fmi2mspct	<i>Read multiple solar spectra from a data file.</i>
----------------	--

Description

Read spectral irradiance file as output by Anders Lindors' model based on libRadTrans for hourly simulation, or measured data from FMI's Brewer spectrometer.

Usage

```
read_fmi2mspct(
  file,
  scale.factor = 0.001,
  geocode = NULL,
  what.measured = NULL,
  how.measured = NULL,
  date.field = 2L,
  time.field = 3L,
```

```

date.format = "ymd",
time.format = "hms",
tz = NULL,
time.shift.min = 0,
locale = readr::default_locale(),
.skip = 0,
.n_max = -1
)

```

Arguments

file	Either a path to a file, a connection, or literal data (either a single string or a raw vector).
scale.factor	numeric A multiplier to be applied to the spectral irradiance values.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
what.measured	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
how.measured	character string, but if NULL or NA the "how.measured" attribute is not set.
date.field, time.field	integer. Word positions in the header line.
date.format	character string. One of "ymd", "ydm", "dmy", or "mdy".
time.format	character string. One of "hms", "hm".
tz	character Time zone used for interpreting times saved in the file header.
time.shift.min	numeric. Time shift with respect to TZ in minutes.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
.skip	Number of lines to skip before reading data.
.n_max	Maximum number of records to read.

Value

read_fmi2mspct() returns a source_mspct object containing source_spct objects as members, time.unit attribute set to "second" and when.measured attribute set to the date-time values extracted from the file body.

Note

See [read_table](#) for details of acceptable values for file. Individual spectra are names based on time and date in ISO format, at the time zone given by tz but the time shift subtracted. Say for times expressed in headers at UTC + 120 min, we use tz = UTC and time.shift.min = 120 to convert times to UTC. This is different from using tz = Europe/Helsinki, which is not invariant through the course of the year because of daylight saving time. Local time zones is not necessarily consistent across years because of changes in legislation. In contrast UTC is more consistent, making it preferable for time series.

read_foreign2mspct	<i>Read multiple foreign files with spectral data</i>
--------------------	---

Description

Read spectra from a homogeneous list of files based on a path and a list of filenames or a path and a search pattern for files. The imported spectra are returned as a single object of one of the collection of spectra classes from package 'photobiology'.

Usage

```
read_foreign2mspct(path = ".", list = NULL, pattern = NULL, .fun, ...)
```

Arguments

path	character A path point to the location of the files.
list	character A vector or list of character strings pointing to files relative to path,
pattern	character A search pattern to select files within path. See list.files which is used internally. Argument ignored is list is non-null.
.fun	function One of the functions exported by this package for reading spectral data.
...	Named arguments passed ot the call to .fun.

Details

This function iterates over a list of file names reading them with the function passed as argument to '.fun' and combines the spectra as a collection of spectra of a class suitable for the spectral objects returned by the argument to '.fun'. This function can either return for each file read either a single spectrum as an object of class 'generic_spct' or a class derived from it, or a collection of spectra of class 'generic_mspct' or a class derived from it. The class of the returned object depends on the class of the member spectra.

Value

An object of class 'generic_mspct' or a class derived from it, containing a collection of member spectra of class 'generic_spct' or of one of the classes derived from it.

read_FReD_csv	<i>Read '.CSV' FReD database.</i>
---------------	-----------------------------------

Description

Reads a CSV data file downloaded from the FReD (Floral Reflectance Database) and then imports wavelengths and spectral reflectance values and flower ID.

Usage

```
read_FReD_csv(
  file,
  date = NA,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)
```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone used for interpreting times saved in the file header.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names. Those relevant should match the format of the CSV file being read.

Value

A reflectance_spct object.

References

<http://www.reflectance.co.uk> Arnold SEJ, Faruq S, Savolainen V, McOwan PW, Chittka L, 2010 FReD: The Floral Reflectance Database - A Web Portal for Analyses of Flower Colour. PLoS ONE 5(12): e14287. doi:10.1371/journal.pone.0014287

Examples

```
file.name <-
  system.file("extdata", "FReDflowerID_157.csv",
             package = "photobiologyInOut", mustWork = TRUE)

fred.spct <- read_FReD_csv(file = file.name)

fred.spct
getWhatMeasured(fred.spct)
cat(comment(fred.spct))
```

read_li180_txt *Read '.TXT' File(s) Saved by LI-COR's LI-180 spectroradiometer.*

Description

Reads and parses the header of a data file as output by the LI-180 spectrometer (not to be confused with the LI-1800 spectrometer released in the 1980's by LI-COR) to extract the whole header remark field and also decode whether data is in photon or energy based units. This is a new instrument released in year 2020.

Usage

```
read_li180_txt(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  s.qty = "s.e.irrad"
)

read_m_li180_txt(
  files,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = Sys.timezone(),
  locale = readr::default_locale(),
  s.qty = NULL
)
```

Arguments

file	Path to file as a character string.
------	-------------------------------------

date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone used for interpreting times saved in the file header.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
s.qty	character The name of the spectral quantity to be read. One of "s.e.irrad" or "s.q.irrad".
files	A list or vector of character strings.

Details

Function `read_m_licor_espd()` calls `red_licor_espd()` for each file in `files`. See [read.table](#) for a description of valid arguments for `files`.

Value

`read_licor_espd()` returns a `source_spct` object with `time.unit` attribute set to "second" and `when.measured` attribute set to the date-time extracted from the file header, or supplied by the user. Spectrometer model, serial number and integration time are stored in attributes. The whole file header is saved as a comment while the footer is discarded.

Function `read_m_licor_espd()` returns a `source_mspct` object containing one spectrum per file read.

Note

The LI-180 spectroradiometer stores little information of the instrument and settings, possibly because they cannot be altered by the user or configured. The length of the file header does not seem to be fixed, so the start of the spectral data is detected by searching for "380nm".

References

LI-COR Biosciences, Environmental. <https://www.licor.com/env/>

Examples

```
file.name <-
  system.file("extdata", "licor-li180-irrad.txt",
             package = "photobiologyInOut", mustWork = TRUE)

licor180.spct <- read_li180_txt(file = file.name)

licor180.spct
getWhenMeasured(licor180.spct)
```

```
getWhatMeasured(licor180.spct)
cat(comment(licor180.spct))
```

read_licor_prn	<i>Read '.PRN' File(s) Saved by LI-COR's PC1800 Program.</i>
----------------	--

Description

Read and parse the header of a processed data file as output by the PC1800 program to extract the whole header remark field and also decode whether data is irradiance in photon or energy based units, transmittance, reflectance or absorbance and then extract the wavelength and spectral data. PC1800 is an MS-DOS program provided for use with the LI-1800 spectrometer. This instrument was released in the 1980's and was sold until the early 2000's. It was very popular and several of them remain in use. (It should not be confused with the LI-180, a new spectrometer released in 2020.)

Usage

```
read_licor_prn(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  s.qty = NULL
)

read_m_licor_prn(
  files,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  s.qty = NULL
)
```

Arguments

file	Path to file as a character string.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.

tz	character Time zone used for interpreting times saved in the file header.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
s.qty	character The name of the spectral quantity to be read. One of "s.irrad", "Tfr", or "Rfr".
files	A list or vector of character strings.

Details

Function `read_m_licor_prn()` calls `read_licor_prn()` for each file in `files`. See [read_table](#) for a description of valid arguments for `files`.

Value

`read_licor_prn()` returns a `source_spct` object with `time.unit` attribute set to "second" and `when.measured` attribute set to the date-time extracted from the file name, or supplied.

Function `read_m_licor_prn()` returns a `source_mspct` object containing one spectrum per file read.

Note

The LI-1800 spectroradiometer does not store the year as part of the data, only month, day, and time of day. Because of this, in the current version, if `NULL` is the argument to `date`, year is set to 0000. In addition, the argument passed to `tz` does not recognize "summer-time" shifts if the year is unknown (date read from the file header).

References

LI-COR Biosciences, Environmental. <https://www.licor.com/env/>

Examples

```
file.name <-
  system.file("extdata", "licor-li1800-q-irrad.PRN",
             package = "photobiologyInOut", mustWork = TRUE)

licor.spct <- read_licor_prn(file = file.name)

licor.spct
getWhenMeasured(licor.spct)
getWhatMeasured(licor.spct)
cat(comment(licor.spct))
```

read_macam_dta	<i>Read '.DTA' File Saved by Macam's Software.</i>
----------------	--

Description

Reads and parses the header of a processed data file as output by the PC program to extract the time and date fields and a user label if present, and then imports wavelengths and spectral energy irradiance values.

Usage

```
read_macam_dta(  
  file,  
  date = NULL,  
  geocode = NULL,  
  label = NULL,  
  tz = NULL,  
  locale = readr::default_locale()  
)
```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone used for interpreting times saved in the file header.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.

Value

A source_spct object.

References

<https://www.irradian.co.uk/>

Examples

```
file.name <-
  system.file("extdata", "macam-e-irrad.DTA",
             package = "photobiologyInOut", mustWork = TRUE)

macam.spct <- read_macam_dta(file = file.name)

macam.spct
getWhenMeasured(macam.spct)
getWhatMeasured(macam.spct)
cat(comment(macam.spct))
```

read_oo_jazirrad

Read Files Saved by Ocean Optics' Jaz spectrometer.

Description

Reads and parses the header of processed data text files output by Jaz instruments extracting the spectral data from the body of the file and the metadata, including time and date of measurement from the header. Jaz modular spectrometers were manufactured by Ocean Optics.

Usage

```
read_oo_jazirrad(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  unit.in = "energy"
)

read_oo_jazpc(
  file,
  qty.in = "Tpc",
  Tfr.type = c("total", "internal"),
  Rfr.type = c("total", "specular"),
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)

read_oo_jazdata(
```

```

    file,
    date = NULL,
    geocode = NULL,
    label = NULL,
    tz = NULL,
    locale = readr::default_locale()
  )

```

Arguments

file	character string Path to the file to be read, following R's use of forward slashes as separator for folder names.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat, and optionally address used to set attribute "where.measured".
label	character string to which to set the "what.measured" attribute. If NULL the value of basename(file) is used, and if NA the "what.measured" attribute is not set.
tz	character A time zone recognized by R. If NULL, the default, it is extracted from 'locale'.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names. Its value must match that used to write the imported file, which is not necessarily the default one or the local one.
unit.in	character One of "energy", "photon" (or "quantum"), for data in uW cm-2 nm-1 and umol cm-2 nm-1.
qty.in	character string, one of "Tpc" (spectral transmittance, %), "A" (spectral absorbance), or "Rpc" (spectral reflectance, %).
Tfr.type	character string, either "total" or "internal".
Rfr.type	character string, either "total" or "specular".

Details

Function read_oo_jazirrad can read processed irradiance output files. Function read_oo_jazpc can read processed transmittance and reflectance output files (expressed as %s). Function read_oo_jazdata can read raw-counts data.

Value

A source_spct object, a filter_spct object, a reflector_spct object or a raw_spct object.

Note

Although the parameter is called date a date time is accepted and expected. Time resolution is < 1 s if seconds are entered with a decimal fraction, such as "2021-10-05 10:10:10.1234".

References

<https://www.oceanoptics.com/>

Examples

```
file.name <-
  system.file("extdata", "oo-spectrum.jaz",
             package = "photobiologyInOut", mustWork = TRUE)

jaz.filter_spct <- read_oo_jazpc(file = file.name)

jaz.filter_spct
getWhenMeasured(jaz.filter_spct)
getWhatMeasured(jaz.filter_spct)
cat(comment(jaz.filter_spct))

file.name <-
  system.file("extdata", "oo-spectrum.jazIrrad",
             package = "photobiologyInOut", mustWork = TRUE)

jaz.source_spct <- read_oo_jazirrad(file = file.name, unit.in = "energy")

jaz.source_spct
getWhenMeasured(jaz.source_spct)
getWhatMeasured(jaz.source_spct)
cat(comment(jaz.source_spct))
q_irrad(jaz.source_spct, waveband(c(400, 700)), scale.factor = 1e6) # mol -> umol
```

read_oo_ovirrad

Read File Saved by Ocean Optics' OceanView.

Description

Reads the spectral data and in addition parses the header of a energy irradiance data file as output by OceanView. OceanView is a program from Ocean Optics used to measure UV, visible and NIR radiation with array spectrometers from the same company. OceanView replaces the no longer supported SpectraSuite program.

Usage

```
read_oo_ovirrad(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  range = NULL
)
```

Arguments

file	character string Path to the file to be read, following R's use of forward slashes as separator for folder names.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat, and optionally address used to set attribute "where.measured".
label	character string to which to set the "what.measured" attribute. If NULL the value of basename(file) is used, and if NA the "what.measured" attribute is not set.
tz	character A time zone recognized by R. If NULL, the default, it is extracted from 'locale'.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names. Its value must match that used to write the imported file, which is not necessarily the default one or the local one.
range	a numeric vector of length two, or any other object for which function range() will return a range of wavelengths expressed in nanometres.

Details

The header of the file is first decoded and parsed to extract the time of data acquisition and serial number of the spectrometer, and to locate the start of the spectral data. The fields are located by name. The spectral irradiance is re-expressed in $Wm^{-2}nm^{-1}$ and returned as an object of class [source_spct](#) with metadata stored in attributes when.measured, what.measured, and how.measured set to values extracted from the header. The value stored in the how.measured attribute includes the User: and Serial Number: values extracted from the file header. If an argument is passed to parameter geocode, its value is saved in attribute where.measured. The file header in whole is copied to attribute file.header. The object's comment always gives a text that includes the file name, time of import, function name and the version of packages 'photobiology' and 'photobiologyInOut' used.

Value

A source_spct object with columns w.length with wavelengths in nanometres and s.e.irrad with spectral energy irradiance in $Wm^{-2}nm^{-1}$, attributes comment, what.measured, when.measured, how.measured, where.measured and file.header containing metadata for the spectrum.

References

<https://www.oceanoptics.com/>

Examples

```
# energy spectral irradiance file
```

```

file.name <-
  system.file("extdata", "oo-spectrum.OV Irrad",
              package = "photobiologyInOut", mustWork = TRUE)

ooov.spct <-
  read_oo_ovirrad(file = file.name,
                  locale = readr::locale("en",
                                          decimal_mark = ",",
                                          grouping_mark = "",
                                          tz = "Europe/Warsaw"))

ooov.spct
getWhenMeasured(ooov.spct)
getWhatMeasured(ooov.spct)
getHowMeasured(ooov.spct)
cat(comment(ooov.spct))

ooov_clipped.spct <-
  read_oo_ovirrad(file = file.name,
                  locale = readr::locale("en",
                                          decimal_mark = ",",
                                          grouping_mark = "",
                                          tz = "Europe/Warsaw"),
                  range = c(280, NA))

ooov_clipped.spct
getWhenMeasured(ooov_clipped.spct)
getWhatMeasured(ooov_clipped.spct)
getHowMeasured(ooov_clipped.spct)
cat(comment(ooov_clipped.spct))

```

read_oo_pidata

Read File Saved by Ocean Optics' Raspberry Pi software.

Description

Reads and parses the header of a raw data file as output by the server running on a Raspberry Pi board to extract the whole header remark field. The time field is retrieved and decoded. The company formerly named Ocean Optics is now called Ocean Insight.

Usage

```

read_oo_pidata(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,

```

```

    locale = readr::default_locale(),
    npixels = Inf,
    spectrometer.sn = "FLMS00673"
  )

```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is set to the file modification date.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone is not saved to the file.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
npixels	integer Number of pixels in spectral data.
spectrometer.sn	character The serial number of the spectrometer needs to be supplied by the user as it is not included in the file header.

Value

A raw_spect object.

Note

The header in these files has very little information. The file contains a time in milliseconds but as the Raspberry Pi board contains no real-time clock, it seems to default to number of milliseconds since the Pi was switched on. The user may wish to supply the date-time as an argument, but if no argument is passed to date this attribute is set to the file modification date obtained with `file.mtime()`. This date-time gives an upper limit to the real time of measurement as in some operating systems it is reset when the file is copied or even without any good apparent reason. The user may need to supply the number of pixels in the array although the default of `npixels = Inf` usually works and triggers no warnings.

References

<https://www.oceanoptics.com/> <https://www.raspberrypi.org/>

Examples

```

file.name <-
  system.file("extdata", "oo-spectrum.pi",
    package = "photobiologyInOut", mustWork = TRUE)

```

```

oopi.spct <- read_oo_pidata(file = file.name)

oopi.spct
getWhenMeasured(oopi.spct)
getWhatMeasured(oopi.spct)
cat(comment(oopi.spct))

```

read_oo_ssirrad

Read File Saved by Ocean Optics' SpectraSuite.

Description

Reads the spectral data and in addition parses the header of a energy irradiance data file as output by SpectraSuite. SpectraSuite is a program from Ocean Optics used to measure UV, visible and NIR radiation with array spectrometers from the same company. OceanView replaces the no longer supported SpectraSuite program.

Usage

```

read_oo_ssirrad(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  range = NULL
)

read_oo_ssdata(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)

```

Arguments

file	character string Path to the file to be read, following R's use of forward slashes as separator for folder names.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat, and optionally address used to set attribute "where.measured".

label	character string to which to set the "what.measured" attribute. If NULL the value of basename(file) is used, and if NA the "what.measured" attribute is not set.
tz	character A time zone recognized by R. If NULL, the default, it is extracted from 'locale'.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names. Its value must match that used to write the imported file, which is not necessarily the default one or the local one.
range	a numeric vector of length two, or any other object for which function range() will return a range of wavelengths expressed in nanometres.

Details

The header of the file is first decoded and parsed to extract the time of data acquisition and serial number of the spectrometer, and to locate the start of the spectral data. The fields are located by name. The spectral irradiance is re-expressed in $Wm^{-2}nm^{-1}$ and returned as an object of class [source_spct](#) with metadata stored in attributes when.measured, what.measured, and how.measured set to values extracted from the header. The value stored in the how.measured attribute includes the User: and Serial Number: values extracted from the file header. If an argument is passed to parameter geocode, its value is saved in attribute where.measured. The file header in whole is copied to attribute file.header. The object's comment always gives a text that includes the file name, time of import, function name and the version of packages 'photobiology' and 'photobiologyInOut' used.

Value

A raw_spct object.

References

<https://www.oceanoptics.com/>

Examples

```
file.name <-
  system.file("extdata", "oo-spectrum.SSIrrad",
             package = "photobiologyInOut", mustWork = TRUE)

ooss.spct <- read_oo_ssirrad(file = file.name)

ooss.spct
getWhenMeasured(ooss.spct)
getWhatMeasured(ooss.spct)
getHowMeasured(ooss.spct)
cat(comment(ooss.spct))
```

read_qtuv_txt	<i>Read Quick TUV output file.</i>
---------------	------------------------------------

Description

Reads and parses the header of a text file output by the Quick TUV on-line web front-end at UCAR to extract the header and spectral data. The time field is converted to a date.

Usage

```
read_qtuv_txt(
  file,
  ozone.du = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  added.vars = NULL
)
```

Arguments

file	character string with the name of a text file.
ozone.du	numeric Ozone column in Dobson units.
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone is by default read from the file.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
added.vars	character vector Accepted member values are "sun.elevation", zenith.angle, "time" and "ozone.du".

Value

a source_spct object obtained by finding the center of wavelength intervals in the Quick TUV output file, and adding the variables listed in added.vars. To obtain the same value as in version (<= 0.4.28) pass added.vars = c("angle", "date") in the call.

Note

The ozone column value used in the simulation cannot be retrieved from the file. Tested files from Quick TUV version 5.2 on 2018-07-30 and also with more recent files in early 2024. This function can be expected to be robust to variations in the position of lines in the imported file and resistant to the presence of extraneous text or even summaries. By default web browsers save the output returned by the Quick TUV calculator as an HTML output, some of them with minimal headers and

other with more extensive ones. In some cases, character escapes replace actual new lines. In most cases these HTML files are decoded correctly, but if not, use "save as" in the browser and select "text" when saving. As a last recourse, messed up files can be manually edited before import.

References

https://www.acom.ucar.edu/Models/TUV/Interactive_TUV/

read_spectrapen_csv *Read '.CSV' File Saved by PSI's Software.*

Description

Reads and parses the header of a processed .CSV file as output by the by the PSI (Photon Systems Instruments, Czech Republic) SpectraPen miniature spectrometer.

Usage

```
read_spectrapen_csv(
  file,
  start.row = 1,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)
```

Arguments

file	character string.
start.row	integer The first line to read, counting from the top of the file.
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured". If NULL, the default, the location is extracted from the file header.
label	character string, but if NULL the label from file header is used, if the label is missing, the index is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone used for interpreting times saved in the file header.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.

Value

A source_mspct object.

References

<https://psi.cz/>

Examples

```
# fetch path to example file to read
file.name <-
  system.file("extdata", "psi-spectrapen-SP.csv",
             package = "photobiologyInOut", mustWork = TRUE)

spectrapen.mspct <- read_spectrapen_csv(file = file.name)

spectrapen.mspct
getWhenMeasured(spectrapen.mspct)
getWhatMeasured(spectrapen.mspct)
cat(comment(spectrapen.mspct))
```

read_tuv_usrout	<i>Read TUV output file.</i>
-----------------	------------------------------

Description

Reads and parses the header of a text file output by the TUV program to extract the header and spectral data. The time field is converted to a date.

Usage

```
read_tuv_usrout(
  file,
  ozone.du = NULL,
  date = lubridate::today(),
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)

read_tuv_usrout2mspct(
  file,
  ozone.du = NULL,
  date = lubridate::today(),
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale()
)
```

Arguments

<code>file</code>	character string
<code>ozone.du</code>	numeric Ozone column in Dobson units.
<code>date</code>	a <code>POSIXct</code> object to use to set the "when.measured" attribute. If <code>NULL</code> , the default, the date is extracted from the file header.
<code>geocode</code>	A data frame with columns <code>lon</code> and <code>lat</code> used to set attribute "where.measured".
<code>label</code>	character string, but if <code>NULL</code> the value of <code>file</code> is used, and if <code>NA</code> the "what.measured" attribute is not set.
<code>tz</code>	character Time zone is by default read from the file.
<code>locale</code>	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.

Value

a `source_spect` object obtained by 'melting' the TUV file, and adding a factor `spect.idx`, and variables `zenith.angle` and `date`.

Note

The ozone column value used in the simulation cannot be retrieved from the file. Tested only with TUV version 5.0.

References

<https://www2.acom.ucar.edu/modeling/tropospheric-ultraviolet-and-visible-tuv-radiation-model>

<code>read_uvspec_disort</code>	<i>Read libRadtran's uvspec output file.</i>
---------------------------------	--

Description

Read and parse a text file output by libRadtran's `uvspec` routine for a solar spectrum simulation. The output of `uvspec` depends among other things on the solver used. We define a family of functions, each function for a different solver.

Usage

```
read_uvspec_disort(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
```

```

    locale = readr::default_locale(),
    multiplier = 0.001,
    qty = "irradiance"
  )

```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone is by default read from the file.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
multiplier	numeric A multiplier for conversion into $\text{W m}^{-2} \text{nm}^{-1}$, as the units of expression of the output from "uvspec" depend on the units in which the extraterrestrial solar spectrum data is expressed.
qty	character "uvspec" returns both irradiance and intensity with solver "disort".

Value

A source_spect object.

Note

Currently only "irradiance" is supported as qty argument as intensity is not supported by classes and methods in package 'photobiology'.

Tested only with libRadtran version 2.0

References

<https://www.r4photobiology.info> <http://www.libradtran.org>

read_uvspec_disort_vesa

Read libRadtran's uvspec output file from batch job.

Description

Reads and parses the header and body of a text file output by a script used to run libRadtran's uvspec in a batch job for a set of solar spectrum simulations. The header and time and date fields are converted into a datetime object.

Usage

```
read_uvspec_disort_vesa(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  multiplier = 1e-06,
  simplify = TRUE
)
```

Arguments

file	character string
date	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date is extracted from the file header.
geocode	A data frame with columns lon and lat used to set attribute "where.measured".
label	character string, but if NULL the value of file is used, and if NA the "what.measured" attribute is not set.
tz	character Time zone is by default read from the file.
locale	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
multiplier	numeric A multiplier for conversion into $\text{W m}^{-2} \text{ nm}^{-1}$, as the units of expression of the output from "uvspec" depend on the units in which the extraterrestrial solar spectrum data is expressed.
simplify	logical Remove redundant columns from returned object.

Value

a source_spect object, possibly containing several spectra in long form and a datetime column.

References

<http://www.libradtran.org>

read_walz_lsa_xlsx	<i>Read '.xlsx' file(s) saved by Walz's LSA-2050 leaf analyser</i>
--------------------	--

Description

Read the two worksheets contained in the workbook file downloaded by Walz's LSA software from an LSA-2050 instrument. Combine the data into a single data frame and add metadata as its attributes.

Usage

```
read_walz_lsa_xlsx(
  file,
  label = NULL,
  data_skip = 0,
  n_max = Inf,
  na = c("", "NA", "NAN"),
  marker.rename = NULL,
  returned.data = "consolidated",
  drop.positions = NULL,
  ...
)
```

Arguments

<code>file</code>	Path to file as a character string.
<code>label</code>	character string, but if <code>NULL</code> the value of <code>file</code> is used, and if <code>NA</code> the "what.measured" attribute is not set.
<code>data_skip</code>	integer Number of records (rows) to skip from the actual data block.
<code>n_max</code>	integer Maximum number of records to read.
<code>na</code>	character Vector of strings to interpret as missing values. Set this option to <code>character()</code> to indicate no missing values.
<code>marker.rename</code>	character A named vector with new Marker values named according to the Markers set in the LSA-2050.
<code>returned.data</code>	character If "consolidated" return all data in a single data frame, with the "SAT chart" data as a list column containing one data frame per sample. Other values return a data frame containing a subset of the columns: with "measured", return only the data in the worksheet named "Measured", with "visible" return only data from columns visible in the worksheet named "Measured", and "raw" returns only data stored in both worksheets as numeric values, i.e., ignores values computed by formulas included in the worksheet. <i>For debugging:</i> if "list", return a named list with two member data frames, each containing the data from one worksheet.
<code>drop.positions</code>	logical If TRUE the columns related to leaf position are deleted, and if FALSE they are retained. If NULL they are dropped only if they contain no finite data.
<code>...</code>	Further named arguments currently passed to <code>read_excel()</code> .

Details

The worksheets contain formulas for computations, on import the results of the computations are computed and stored in the returned object. The raw readings and calibration values are also imported as stored as an attribute. The names of the worksheets are preserved, and if they do not match the expected values, a warning is issued. Data manipulation is based on column names rather than positions. As the Date and Time columns contain identical data (formatted differently) for efficiency, only Time is retained in the returned object (after testing for equality). Units and bases of expression are not modified, but some column names are edited to comply with R naming rules.

The named vector passed as argument to `marker.rename` if not `NULL`, is used to rename the values set as Marker by the LSA-2050 (single capital letters) into informative character strings. `marker.rename` can contain a superset of the required mappings but a subset is not accepted.

Value

`read_csi_dat()` returns a `tibble::tibble` object.

The returned value depends on the argument passed to parameter `returned.data`. With the default, `"consolidated"`, all the data from both worksheets are returned in a single data frame with 44 columns. With arguments `"measured"` (43 columns), `"visible"` (17 columns) and `"raw"` (24 columns) different variable subsets from the consolidated data frame are returned. Finally for debugging, with `"list"`, a list containing one data frame with the data from each worksheet is returned.

Warning!

The computed values for the formulas in the worksheets need to be present in the Excel workbook for import to succeed. If they are missing, zeros will appear in the worksheet in place of the results from calculations. The workbooks saved using the LSA software contain the computed cached values. Opening and re-saving the file can result in the loss of the cached computed values. Importing such a file triggers a message. The cached values from formulas can be restored in Excel or in LibreOffice. In LibreOffice open the workbook file, use **data > calculate > recalculate hard** to recalculate all formulas and then save the file. (In LibreOffice a manually triggered hard recalculation is also needed for the correct computed values to be displayed.) The example file used in the examples and tests has been resaved in LibreOffice.

To avoid depending on the import of cached values stored in the worksheets, one possible solution would be to reimplement all computations in R after importing the `"raw"` data (plus calibration parameters). These are stored in the worksheet as numeric values. This would provide a failure-proof approach, but updates to equations used in the worksheets would require matching updates to the R code. Passing `returned.data = "raw"` in the call imports only these values.

Note

This function has been tested with the output from LSA software version 1.09 and an example data file downloaded from Walz's website.

References

<https://www.campbellsci.eu/>

Examples

```
file.name <-
  system.file("extdata", "walz-lsa-2050-prunus.xlsx",
             package = "photobiologyInOut",
             mustWork = TRUE)

lsa.df1 <- read_walz_lsa_xlsx(file.name)
colnames(lsa.df1)
plot(lsa.df1$SAT.F.ls[[1]], type = "b")
```

```

cat(comment(lsa.df1))

lsa.df2 <-
  read_walz_lsa_xlsx(file.name,
                     returned.data = "visible",
                     marker.rename = c(A = "adaxial", B = "abaxial"))
colnames(lsa.df2)
cat(comment(lsa.df2))

```

read_wasatch_csv

Read File Saved by Wasatch's Enlighten.

Description

Read wavelength and spectral data from the data section of a file as output by Enlighten importing them into R. Parse the header of a file to extract the acquisition time, instrument name and serial number, as well additional metadata related to the instrument and its settings. Function `read_wasatch_csv()` only accepts "column oriented" CSV files.

Usage

```

read_wasatch_csv(
  file,
  date = NULL,
  geocode = NULL,
  label = NULL,
  tz = NULL,
  locale = readr::default_locale(),
  s.qty = NULL,
  extra.cols = "keep",
  scale.factor = 1,
  simplify = TRUE,
  ...
)

```

Arguments

<code>file</code>	character
<code>date</code>	a POSIXct object to use to set the "when.measured" attribute. If NULL, the default, the date and time are extracted from the file header.
<code>geocode</code>	A data frame with columns <code>lon</code> and <code>lat</code> used to set attribute "where.measured".
<code>label</code>	character string, but if NULL the value of <code>file</code> is used, and if NA the "what.measured" attribute is not set.
<code>tz</code>	character Time zone is by default that of the machine's locale.

<code>locale</code>	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use <code>locale</code> to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.
<code>s.qty</code>	character, possibly named. The name of the quantity using the conventions accepted used in package 'photobiology' that is to be imported from column "Processed" from the file.
<code>extra.cols</code>	character What to do non-processed data columns if present in file. One of "keep", "drop.pixel", "drop" or "split".
<code>scale.factor</code>	numeric vector of length 1, or length equal to the number of rows (= detector pixels). Numeric multiplier applied to returned spectral values.
<code>simplify</code>	logical If TRUE, single spectra are returned as individual spectra instead of collections of length one.
<code>...</code>	additional arguments passed to the constructor of the spectrum object.

Details

Enlighten's column-wise CSV files contain at least two columns, Wavelength and Processed. In the header the Technique used is recorded. Additional data columns can be present. Column Pixel contains the pixel index in the array as integers. Columns Raw, Dark and Reference contain detector counts data. Technique is used to guess the type of spectrum stored in the column named Processed, which can be detector counts or derived values. By default the data are read into a single spectrum object and all columns retained, but only the data in Processed are interpreted as spectral data corresponding to the class of the object. If passed `extra.cols = "drop"`, only Wavelength and Processed are copied to the returned object, while if passed `extra.cols = "drop.pixel"` only the contents of column Pixel are discarded. If passed `extra.cols = "split"` all columns containing spectral data are each read into a separate spectrum, these are collected and a "generic_mspect" object containing them returned. `extra.cols` can be a named vector of mappings, of length at least one but possibly longer. If longer a "generic_mspect" is returned, otherwise a spectrum object as inferred from the name each column is mapped to.

Value

An object of a class derived from `generic_spct` such as `raw_spct` or `filter_spct`. `generic_spct` is derived from tibble and data frame.

Acknowledgements

We thank Ruud Niesen from Photon Mission (<https://photonmission.com/>) for organizing the loan of the spectrometer used to produce the various files needed for the development of this function.

Note

Enlighten, the free software from Wasatch Photonics can save spectra in a variety of additional formats: different types of CSV files, plain text and JSON. Plain text files contain no metadata or even column headers and if the need arises can be read with R function `read.table()`. JSON files contain the most detailed metadata.

References

<https://wasatchphotonics.com/> <https://wasatchphotonics.com/product-category/software/>

Examples

```
file.name <-
  system.file("extdata", "wasatch-enlighten-scope.csv",
             package = "photobiologyInOut", mustWork = TRUE)

wasatch.raw.spct <-
  read_wasatch_csv(file = file.name)
summary(wasatch.raw.spct)

wasatch.raw.spct <-
  read_wasatch_csv(file = file.name, s.qty = "counts")
summary(wasatch.raw.spct)

wasatch.raw.spct <-
  read_wasatch_csv(file = file.name, s.qty = c(Processed = "counts"))
summary(wasatch.raw.spct)

wasatch.raw.spct <-
  read_wasatch_csv(file = file.name, extra.cols = "drop")
summary(wasatch.raw.spct)
```

read_yoctopuce_csv	<i>Read '.CSV' file(s) downloaded from YoctoPuce modules.</i>
--------------------	---

Description

Reads and parses the header of processed data CSV files as output by the virtual- or hardware-hubs and modules from Yoctopuce. Uses the comment attribute to store the metadata.

Usage

```
read_yoctopuce_csv(
  file,
  geocode = NULL,
  label = NULL,
  data_skip = 0,
  n_max = Inf,
  locale = readr::default_locale()
)
```

Arguments

<code>file</code>	Path to file as a character string.
<code>geocode</code>	A data frame with columns <code>lon</code> and <code>lat</code> used to set attribute <code>"where.measured"</code> .
<code>label</code>	character string, but if <code>NULL</code> the value of <code>file</code> is used, and if <code>NA</code> the <code>"what.measured"</code> attribute is not set.
<code>data_skip</code>	integer Number of records (rows) to skip from the actual data block.
<code>n_max</code>	integer Maximum number of records to read.
<code>locale</code>	The locale controls defaults that vary from place to place. The default locale is US-centric (like R), but you can use locale to create your own locale that controls things like the default time zone, encoding, decimal mark, big mark, and day/month names.

Details

Yoctopuce modules are small USB connected and USB powered, but isolated, very high quality miniature data acquisition and interface modules. All modules capable of data acquisition can log measured data autonomously and these data can be locally or remotely downloaded as a CSV file. (It is also possible and very easy to access these modules from R using package 'reticulate' and the Python library provided by Yoctopuce, or to send commands and retrieve data through the built-in HTML server of the modules or dedicated hubs.)

Value

`read_yoctopuce_csv()` returns a `tibble::tibble` object, with the number of columns dependent on the CSV file read.

Note

This function should be able to read data log files from any YoctoPuce USB interface module with data logging capabilities as the format is consistent among them.

References

<https://www.yoctopuce.com/>

Examples

```
# We read a CSV file previously downloaded from a YoctoMeteo module.

file.name <-
  system.file("extdata", "yoctopuce-data.csv",
             package = "photobiologyInOut", mustWork = TRUE)

yoctopc.tb <- read_yoctopuce_csv(file = file.name)

yoctopc.tb
cat(comment(yoctopc.tb))
```

rspec2mspct	<i>Convert "pavo::rspec" objects</i>
-------------	--------------------------------------

Description

Convert between objects of class `rspec` as used in package 'pavo' and similar objects returned by functions from package 'lightr' in objects of classes from 'photobiology'

Usage

```
rspec2mspct(
  x,
  member.class = "reflector_spct",
  spct.data.var = "Rpc",
  multiplier = 1,
  ...
)

rspec2spct(x, multiplier = 1, ...)
```

Arguments

<code>x</code>	<code>rspec</code> object
<code>member.class</code>	character One of the spectrum classes defined in package 'photobiology'. Examples are 'reflector_spct', 'filter_spct', 'source_spct', etc.
<code>spct.data.var</code>	character The name to be used for the 'spc' data when constructing the spectral objects. Examples are spectral reflectance as percent or as a fraction of one, 'Rpc', 'Rfr', spectral transmittance 'Tpc', 'Tfr', spectral energy irradiance in Wm^{-2} 's.e.irrad', spectral photon irradiance in $molm^{-2}s^{-1}$'s.q.irrad', etc.
<code>multiplier</code>	numeric A multiplier to be applied to spectral data in the 'rspec' object to do unit or scale conversion.
<code>...</code>	currently ignored.

Details

`rspec` objects contain data for multiple spectra in wide form, i.e., one column per spectrum. Function `rspec2mspct()` converts them into a list-like collection of spectra (`xxxx_mspct`) with spectra stored as separate members. In contrast, `rspec2spct()` stores the spectra in long form into a single data-frame-like object (`xxxx_spct`). The classes of the returned objects are defined in package 'photobiology' and derived from `list` and `data.frame` respectively.

Objects of class `pavo::rspec` do not contain metadata or class data from which the quantity measured and the units of expression could be retrieved. Thus, when using this function users need to use parameter `multiplier` to convert the data to the units and scaling that is expected by the object constructors defined in package 'photobiology' and use parameter `spct.data.var` to indicate the physical quantity.

@note `pavo::rspec` objects may use memory more efficiently than spectral objects of the classes for collections of spectra defined in package 'photobiology' as in 'pavo' wavelengths are assumed to be the same for all spectra stored in an object. So, wavelengths are always stored only once per object. In contrast, 'photobiology' avoids this assumption allowing heterogeneous collections containing member spectra with different wavelengths and lengths. This requires storing Wavelengths for each spectrum, but as spectral classes are derived from 'tbl_df' in many cases redundant copies of wavelength data are avoided by optimizations in 'tibble'.

Warning!

Always check the sanity of the imported or exported data values, as guessing is needed when matching the different classes, and the functions defined here are NOT guaranteed to return valid data without help from the user through optional function arguments.

Examples

```
# example run only if 'pavo' is available
if (requireNamespace("pavo", quietly = TRUE)) {
  library(pavo)
  data(sicalis, package = "pavo")
  sicalis.mspct <- rspec2mspct(sicalis)
  class(sicalis.mspct)

  data(teal, package = "pavo")
  teal.spct <- rspec2spct(teal)
  class(teal.spct)
  levels(teal.spct[["spct.idx"]])
  angles <- seq(from = 15, to = 75, by = 5) # from teal's documentation
  teal.spct[["angle"]] <- angles[as.numeric(teal.spct[["spct.idx"]])]
  teal.spct
}
```

spct_CCT

Correlated color temperature

Description

Wrapper on function `computeCCT` from package 'colorSpec' that accepts `source_spct` objects.

Usage

```
spct_CCT(
  spct,
  isotherms = "robertson",
  locus = "robertson",
  strict = FALSE,
  named = FALSE
)
```

Arguments

spct	source_spct A single light source spectrum.
isotherms	character A vector whose elements match one of the available isotherm families: 'robertson', 'mccamy', and 'native'. Matching is partial and case-insensitive. When more than one family is given, a matrix is returned, see Value. When isotherms = 'native' the isotherms are defined implicitly as lines perpendicular to the locus, see Details in CCTfromXYZ. The character NA (NA_character_) is taken as a synonym for 'native'.
locus	character Valid values are 'robertson' and 'precision', see above. Matching is partial and case-insensitive.
strict	logical The CIE considers the CCT of a chromaticity uv to be meaningful only if the distance from uv to the Planckian locus is less than or equal to 0.05 (in CIE UCS 1960). If strict=FALSE, then this condition is ignored. Otherwise, the distance is computed along the corresponding isotherm, and if it exceeds 0.05 the returned CCT is set to NA.
named	logical Whether to set the name attribute of the returned value to the name of the spectrum passed as argument if possible.

Details

Please see [computeCCT](#) for the details of the computations and references.

Value

A numeric value for "color temperature " in degrees Kelvin.

Examples

```
spct_CCT(white_led.source_spct)
spct_CCT(sun.spct)
```

spct_CRI	<i>Color reproduction index</i>
----------	---------------------------------

Description

Wrapper on function [computeCRI](#) from package 'colorSpec' that accepts [source_spct](#) objects.

Usage

```
spct_CRI(spct, adapt = TRUE, attach = FALSE, tol = 0.0054, named = FALSE)
```

Arguments

spct	source_spct A single light source spectrum.
adapt	logical If TRUE, then a special chromatic adaption is performed, see Details in computeCRI .
attach	logical If TRUE, then a large list of intermediate calculations is attached to the returned number, as attribute data. This attached list includes data for all special 14 color samples, although the last 6 do not affect the returned CRI
tol	numeric For the CRI to be meaningful the chromaticities of the test and reference illuminants must be sufficiently close in the CIE 1960 uniform chromaticity space. If the tolerance is exceeded, the function returns NA. The default tol=5.4e-3 is the one recommended by the CIE, but the argument allows the user to override it.
named	logical Whether to set the name attribute of the returned value to the name of the spectrum passed as argument if possible.

Details

Please see [computeCRI](#) for the details of the computations and references.

Value

A numeric value between zero and 100, or NA if the light is not white enough.

Examples

```
spct_CRI(white_led.source_spct)
spct_CRI(sun.spct)
```

spct_SSI	<i>Spectral (color) similarity index</i>
----------	--

Description

Wrapper on function [computeSSI](#) from package 'colorSpec' that accepts [source_spct](#) objects.

Usage

```
spct_SSI(
  spct,
  reference.spct = NULL,
  digits = 0,
  isotherms = "mccamy",
  locus = "robertson",
  named = FALSE
)
```

Arguments

spct, reference.spct	source_spct Single light source spectra.
digits	integer The number of digits after the decimal point in the returned vector. According to Holm the output should be rounded to the nearest integer, which corresponds to <code>digits = 0</code> . To return full precision, set <code>digits = Inf</code> .
isotherms	character This is only used when <code>reference=NULL</code> . It is passed to computeCCT in order to compute the CCT of each test spectrum.
locus	character This is only used when <code>reference=NULL</code> . It is passed to computeCCT in order to compute the CCT of each test spectrum.
named	logical Whether to set the name attribute of the returned value to the name of the spectrum passed as argument if possible.

Details

Please see [computeSSI](#) for the details of the computations and references.

Value

A numeric value between zero and 100.

Examples

```
spct_SSI(white_led.source_spct, sun.spct)
```

Index

* misc

- read_avaspec_csv, 19
- read_CIE_csv, 23
- read_FReD_csv, 30
- read_licor_prn, 33
- read_oo_ovirrad, 38

- as.chroma_mspct.colorSpec
 (colorSpec2mspct), 4
- as.chroma_spct.colorSpec
 (colorSpec2mspct), 4
- as.colorSpec (mspct2colorSpec), 7
- as.filter_mspct.colorSpec
 (colorSpec2mspct), 4
- as.filter_spct.colorSpec
 (colorSpec2mspct), 4
- as.generic_mspct.colorSpec
 (colorSpec2mspct), 4
- as.generic_spct.colorSpec
 (colorSpec2mspct), 4
- as.reflector_mspct.colorSpec
 (colorSpec2mspct), 4
- as.reflector_spct.colorSpec
 (colorSpec2mspct), 4
- as.response_mspct.colorSpec
 (colorSpec2mspct), 4
- as.response_spct.colorSpec
 (colorSpec2mspct), 4
- as.source_mspct.colorSpec
 (colorSpec2mspct), 4
- as.source_spct.colorSpec
 (colorSpec2mspct), 4

- CCTfromXYZ, 58
- chroma_spct2colorSpec
 (mspct2colorSpec), 7
- colorSpec2chroma_spct
 (colorSpec2mspct), 4
- colorSpec2mspct, 4
- colorSpec2spct (colorSpec2mspct), 4

- computeCCT, 57, 58, 60
- computeCRI, 58, 59
- computeSSI, 59, 60

- fdata2mspct (mspct2fdata), 9
- fdata2spct (mspct2fdata), 9
- filter_spct, 15

- generic_spct, 15

- hyperSpec2mspct, 6
- hyperSpec2spct (hyperSpec2mspct), 6

- list.files, 29
- locale, 12, 15, 18, 20, 21, 25, 26, 28, 30, 32,
 34, 35, 37, 39, 41, 43–45, 47–49, 53,
 55

- mspct2colorSpec, 7
- mspct2fdata, 9
- mspct2hyperSpec (hyperSpec2mspct), 6

- photobiologyInOut
 (photobiologyInOut-package), 3
- photobiologyInOut-package, 3

- qtuv_clouds, 10
- qtuv_m_s.e.irrad (qtuv_s.e.irrad), 11
- qtuv_s.e.irrad, 10, 11

- raw_spct, 15
- read.table, 32
- read_asd_txt, 14
- read_ASTER_txt, 17
- read_avaspec_csv, 19
- read_avaspec_xls (read_avaspec_csv), 19
- read_cid_spectravue_csv, 20
- read_CIE_csv, 23
- read_csi_dat, 24
- read_excel, 50
- read_fmi2mspct, 27

read_fmi_cum, [26](#)
read_foreign2mspct, [29](#)
read_FReD_csv, [30](#)
read_li180_txt, [31](#)
read_licor_prn, [33](#)
read_m_fmi_cum(read_fmi_cum), [26](#)
read_m_li180_txt(read_li180_txt), [31](#)
read_m_licor_prn(read_licor_prn), [33](#)
read_macam_dta, [35](#)
read_oo_jazdata(read_oo_jazirrad), [36](#)
read_oo_jazirrad, [36](#)
read_oo_jazpc(read_oo_jazirrad), [36](#)
read_oo_ovirrad, [38](#)
read_oo_pidata, [40](#)
read_oo_ssdata(read_oo_ssirrad), [42](#)
read_oo_ssirrad, [42](#)
read_qtuv_txt, [44](#)
read_spectrapen_csv, [45](#)
read_table, [27](#), [28](#), [34](#)
read_tuv_usrout, [46](#)
read_tuv_usrout2mspct
 (read_tuv_usrout), [46](#)
read_uvspec_disort, [47](#)
read_uvspec_disort_vesa, [48](#)
read_walz_lsa_xlsx, [49](#)
read_wasatch_csv, [52](#)
read_yoctopuce_csv, [54](#)
reflector_spct, [15](#)
rspec2mspct, [56](#)
rspec2spct(rspec2mspct), [56](#)

source_spct, [15](#), [39](#), [43](#), [57–59](#)
spct2colorSpec(mspct2colorSpec), [7](#)
spct2fdata(mspct2fdata), [9](#)
spct2hyperSpec(hyperSpec2mspct), [6](#)
spct_CCT, [57](#)
spct_CRI, [58](#)
spct_SSI, [59](#)